




# GS MEDICAL PACKAGING INC.

30 Years of Advanced Solutions in all of your Infection Control Needs




All STERIL-PEEL\* product are manufactured from medical grade bleached kraft paper. Hydrogen peroxide, rather than chlorine is used as the bleaching agent.


They are made from virgin pulp, ensuring optimum consistency of the paper sheet, providing superior barrier properties.

STERIL-PEEL\* self-sealing pouches are manufactured with a pre-fold (see diagram) to ensure ease and accuracy of pouch sealing.

The adhesive tape is made of a polyester film coated with medical grade acrylic adhesive, designed for use in steam and ethylene oxide sterilizers.

## QUALITY MANUFACTURING

- 

Compatible with all of the most commonly used sterilization methods. Validated by independent laboratories to ensure that all parameters required have been met in compliance with ISO standards.
- 

Easy identification of contents through transparent blue film. Blue film vs white film – Blue film being a stronger identifier of defects, punctures and/or tears in film.

Water-based adhesives  
Lead-Free Inks  
Dual Indicators  
Latex-free

GS Medical Packaging enjoys over 30 years of being the pioneer in high quality medical pouches and sterilization products. We are one of the few companies that manufacture self-seal pouches with Class 4 indicators. Since our inception, we have created packages that are a cut above the rest in strength, sterility and protection for medical devices, components, trays and cassettes.

Our superiority in design, quality and function is indisputable. These are trademarks of our brand.

GS Medical Packaging's reputation continues to grow both Nationally and Internationally due to our exceptional products and our extensive in-house knowledge of infection control.